

# **Remote User Guide**





This document is provided as a supplement to the Instructions for Use. Always refer to the Instructions For Use for complete operating instructions, warnings and precautions. Avail Medsystems, Inc. Santa Clara, CA 95054 avail.io



# **Table of Contents**

Overview	
Avail Console	4,5
Remote App	
Quick Reference	7
Prepare for Case	8
Remote User Setup	
Pre-Call Setup	
Remote App Setup	9
Log On	9
Accepta Call	9
Set Up Interface	
Collaboration Features	
Views	
Camera Controls	11
Presents	
Freeze / Unfreeze	11
Refresh	11
Annotate / Erase	11
Multiparty Events	12
Broadcasting a Case	13
Troubleshooting	15

#### APPENDIX

Remote User   Initial Setup	. 19
Remote User   Avail Portal	20
Multiparty Events   Setup	. 22
Multiparty Events   Joining	23

#### **Initial Setup**

First time Remote Users will need to complete the initial setup section found in the appendix.



# **Overview**

The Avail System and services are intended to passively transfer and display information for remote sharing and are not intended to diagnose, treat, cure, or prevent any disease.

### **REMOTE USER**

REMOTE APP Enables on-the-go collaboration with remote network members



#### AVAIL PORTAL Allows Remote Users to schedule availability by facility





#### **PROCEDURE ROOM**

AVAIL CONSOLE Provides access to medical procedures anytime, anywhere



AVAIL PORTAL Allows facility admins to control Remote User access

# **Technical Support**



#### 1-833-GO-AVAIL (462-8245) Available 24/7/365

The Avail Support team is available 24 hours a day, seven days a week for facilities if the Avail Console is not running as expected.

#### Internet Bandwidth (WAN)

Test prior to starting a call to ensure consistent media quality. Poor internet speeds may result in reduced media quality.

#### **Recommended Bandwidth**

*at Remote User Location* 20Mbps; minimum speed of 5Mbps



# Avail Console, AC-150





# Avail Console, AC-200





# **Remote App**

### **Overview**





# **Quick Reference**

### **Prepare for Case**



### POINTS OF CONTACT Determine support team

### **Remote User Setup**



### PRE CALL SETUP

Conduct speed test and prepare device (fully charged and hardwire ethernet connection)

### **Remote App Setup**

LOG ON Use either iPad or laptop

### ACCEPT A CALL

Tap "Answer"; View upcoming appointments in calendar

#### **SET UP INTERFACE**

Set picture-in-picture, camera and audio options

### **Collaboration Features**

### VIEWS

Decide display view(s) and confirm external inputs are working

#### **CAMERA CONTROLS**

Pan/zoom front and overhead cameras

#### PRESETS

Save preset; presets are only saved for that Console and remote device

#### **FREEZE/UNFREEZE**

Tap the icon to freeze/resume the screen

#### REFRESH

Tap the icon to refresh the screen

#### Noise Reduction

J.

#### **NOISE REDUCTION**

Remove distracting OR background noise

#### ANNOTATE/ERASE Highlight or draw on s

Highlight or draw on screen; erase as needed

#### **END CALL**

Tap "End Call" icon and close browser



# **Prepare for Case**



#### **POINTS OF CONTACT**

- Determine Avail support point person
- Determine Facility point person

# **Remote User Setup**



#### PRE-CALL SETUP

#### Speed Test

- Using approved browser go to <a href="http://networktest.avail.io/">http://networktest.avail.io/</a> to ensure 20Mbps is met for both upload & download speed
- · Ensure other devices on network are not consuming bandwidth

#### **Prepare Device**

- · Confirm device is fully charged and connected to AC power
- · Use hardwire Ethernet connection & turn off WiFi



#### iPad

- Close all apps
- Turn off VPN

#### **Device Compatibility:**

- iPad OS v12.3 or above
- iPad: 4th, 5th, 6th & 7th Gen
- iPad Air, Air2 & Air3
- iPad Pro: 2nd, 3rd & 4th Gen
- iPad Mini is not currently supported

Atom	
Lag b Republication Ra	
No.	
Andrew Constanting Systems	

#### Laptop

- Close all programs except for compatible browser
- Clear browser cache
- Turn off VPN

#### System Requirements:

- Mac: Mac OS 11+
- Windows: Windows 10
- Microsoft Surface Pro: Windows 10S
   or 10 Pro

#### **Compatible Browsers:**

- Google Chrome (best practice)
- Microsoft Edge
- · Safari (Mac only)
- Firefox



# **Remote App Setup**



#### LOG ON

• Use either iPad or Laptop to access

Open Remote App

Tap "login"

• The same login credentials are used for both

Enter username and password

NOTE If inactive for > 30 minutes, the system will automatically logout.



Go to avail.io Click "Member Login" upper right corner Enter username and password



#### ACCEPT A CALL

- View upcoming appointments & missed calls in calendar
- If a call is missed, you have 15 minutes to call the Facility back on that Console



Tap "Answer" on the incoming call from the Remote App



Click "Answer" on the incoming call from the Portal *A new browser tab will open* 

**TIP:** Monitor internet by periodically confirming five green bars for the internet connection throughout the case

# SET UP USER INTERFACE

**TIP**: If you are Broadcasting, see page 13





Tap Camera icon to turn camera on/off



Tap Microphone icon to mute or unmute

#### To Maximize Browser

Windows: Select F11 Mac: Command, Control, F

#### Chrome Users:

Go to "View", confirm these are NOT selected:

- "Always Show Bookmarks Bar" "Always Show Toolbar in Full
- Screen"
- "Always Show Full URLs"



# **Collaborative Features**

#### 



### VIEWS

- · Decide on the correct display view(s) for procedure
  - Tap the bottom left corner display icons to set which views to display
- · Confirm external inputs are working properly
  - Place in split screen display, if desired



### **CAMERA CONTROLS**

only for front & overhead camera



Zoom In: Double-tap anywhere on the screen to move camera, and use the slider bar to zoom in & out *Using the front camera in full screen* 



Zoom In: Click camera control icon to pan and tilt, and use the slider bar to zoom in & out

# +++

### PRESETS

- Tap the preset icon and select save new preset
- · Enter desired name and tap save



### FREEZE/UNFREEZE

- Tap icon to freeze screen
- Tap the same icon to resume

NOTE

Presets are only saved for that Console and remote device and will need to be reset at the beginning of each call.



#### REFRESH

- Tap icon to refresh the screen
- DO NOT use browser refresh button

#### ANNOTATE/ERASE

- Tap on any of the four colored circles to choose color (contrasting colors work best)
- Highlight or draw on screen using mouse, fingertip, or Apple Pencil
- Tap eraser icon to clear all annotations



#### END CALL

- Tap "End Call" icon
- Close browser



# **Multiparty Events**

A multiparty event allows a Console and up to four Remote Users to connect and collaborate at the same time.



### HOST A MULTIPARTY EVENT

#### Create Event

- Login to Avail Portal at avail.io
- □ Click "Events"
- Click "Create Event"
- □ Enter information and submit

#### Start Event

Event must be initiated by Console

- □ From main Console page, select "Events"
- Select the Event you wish to start

#### **Collaborate During Event**

- Host controls views, camera, annotations, sidebar and audio/video muting
- Use Sidebar mode to converse with the Remote Users only, without disrupting the Console user



#### JOINING AS A REMOTE USER

#### Join via Member Portal

- □ Login to Avail Portal at <u>avail.io</u>
- □ Click "Events"
- □ Find the Event and click "Join"

#### **Tips & Tricks**

- Events can be scheduled in advance and can be joined at any time on the day of the event.
- Events can be edited after creation.
- Non-member guests can be invited to an event and can join as participants.
- During the event, ensure all other applications that use the camera are closed.

#### Join via Event Notification Email

- Open email
- □ Click "Join Now"



# **Broadcasting a Case**

A live broadcast event can occur when a Remote Moderator shares their laptop screen with a live audience using any third-party HIPAA compliant video conferencing platform.

To host a broadcasting event, follow these instructions in addition to all prior steps:



#### **IDENTIFY SUPPORT**

- Determine Broadcast Team
- Assign co-host to monitor the chat for questions
- □ Have a backup moderator ready and available



#### PREPARE VIDEO CONFERENCING PLATFORM

# Confirm General Settings

 $\hfill\square$  Disable local and cloud recording capabilities

#### **Create Meeting/Webinar**

#### Suggested Settings

- □ Mute participants upon entry
- Required Registration
- □ Enable Waiting Room

#### **Create Meeting Invite**

Copy and paste registration link and send to attendees

#### Schedule Test Call (day before or day of)

□ Conduct test call with facility duplicating the meeting requirements for live call



# Broadcasting a Case (continued)



#### **PRE-BROADCAST SETUP**

#### **Network Speed Test for Broadcasting**

□ Using approved browser go to <u>http://networktest.avail.io/</u> to view download speeds (Minimum 100Mbps download and 20Mbps upload speeds)

#### **Prepare Laptop**

- Confirm laptop & browser meet Avail requirements
- Use hard-wired ethernet connection & turn off Wi-Fi
- Close all programs except for compatible browser and video conferencing software
- □ Turn off VPN
- Clear browser cache
- Turn on Do Not Disturb
- Accept Avail Console call before connecting to video conferencing software



#### START VIDEO CONFERENCING CALL

#### Share Screen

- □ Maximize browser
- □ Select "share sound" and "optimize for video clip"
- □ Select window for sharing

#### **Interface Setup**

Set up audio and video per the table below

#### **Conduct Broadcast**

Accept waiting room attendees

#### **To Maximize Browser**

Windows: Select F11 Mac: Command, Control, F

#### **Chrome Users:**

Go to "View", confirm these are NOT selected:

- "Always Show Bookmarks Bar"
- "Always Show Toolbar in Full Screen"
- "Always Show Full URLs"

	<u><u><u></u></u></u>	Q		0
Co	onsole Re Audio Use	emote er Audio P	Picture in Picture (PiP)	Remote User Camera
Broadcast N	IUTE N	<b>/</b> UTE	OFF	OFF
Broadcast M	Audio Use	AUTE	Picture (PiP) OFF	Camera



# Troubleshooting

### **REMOTE APP: iPad CALL**

#### **User account locked**

ŝ	User entered incorrect credentials 5 times in a row	no	Attempt login after 15 min or reset password by selecting "Forgot Password" on the login screen
Cai	Temporary password incorrectly entered 3 times (or) account locked for security reasons	Act	Contact Avail Technical Support

#### Login failed

ause	Potential issue with Remote User's internet connectivity	Action	Check internet speed using a supported browser by going to <u>http://networktest.avail.io/</u> Recommended internet speed is 20Mbps (minimum 5Mbps)
0	Wrong username or password	A	Check for username and password accuracy

#### **Remote App does not annotate**

Cause	Using Apple pencil with low battery (or) pencil not working because of Bluetooth interference	Action	Use tip of finger to annotate
-------	---	--------	-------------------------------

#### Poor call quality, latency

Cause	Potential issue with Remote User's internet connectivity	Action	Check internet speed using a supported browser by going to <u>http://networktest.avail.io/</u> Recommended internet speed is 20Mbps (minimum 5Mbps)
-------	--	--------	--

#### Cannot switch between overhead/front camera or split-screen/full screen views

#### Microphone and/or camera is not working

Cause	iPad does not have permission to access your device's microphone and/or camera	Action	Click "OK" when prompted to grant camera and microphone access to Avail Remote App
-------	--	--------	--



# Troubleshooting

### **REMOTE APP: BROWSER CALL**

#### **User account locked**

ISe	User entered incorrect credentials 5 times in a row	on	Attempt login after 15 min or reset password by selecting "Forgot Password" on the login screen
Cai	Temporary password incorrectly entered 3 times (or) account locked for security reasons	Act	Contact Avail Technical Support

#### Login failed

8	No network connectivity	u	Check network requirements
Cau	Wrong username or password	Acti	Check for username and password accuracy

#### Poor call quality, latency

Cause	Potential issue with Remote User's internet connectivity	Action	Check internet speed using a supported browser by going to <u>http://networktest.avail.io/</u> Recommended internet speed is 20Mbps (minimum 5Mbps)
-------	--	--------	--

#### Cannot switch between overhead/front camera or split-screen/full screen views

Potential software issue	Contact Avail Technical Support
--------------------------	---------------------------------

#### **Firefox browser freezes**

#### **Video appears Pixelated**

Cause	Potential browser rendering issue on certain laptops or desktops	Action	See Help Center for troubleshooting steps by browser type
-------	--	--------	---



# Troubleshooting

### **AVAIL PORTAL**

#### **User account locked**

Cause	User entered incorrect credentials 5 times in a row	Action	Attempt login after 15 min or reset password by selecting "Forgot Password" on the login screen
-------	---	--------	---

#### Login failed

ISe	No network connectivity		Check network requirements
Cau	Wrong username or password	Act	Check for correct username and password entry

#### Remote User cannot schedule calendar appointment for a facility

use	Facility is not listed in remote users' "Facilities"		Select facilities under "My Profile" and click on 'Add Hospital" to request access to the facility
Ca	Facility has not approved Remote User's access request	Ac	Select facilities under "My Profile" and ensure facility access has been approved

#### Scheduled calendar availability appears to be deleted

Cause	Availability likely scheduled as a recurring event and the end date for recurrences has passed	Action	Reschedule availability as recurring event with a new future end date
-------	--	--------	---

#### Scheduled calendar availability appears to be incorrect

Cause	Availability likely scheduled in a different time zone	Action	Check time zone displayed in the drop- down list on the top right; ensure correct time zone is selected
-------	--	--------	---



# **APPENDIX**

Remote User   Initial Setup	19
Remote User   Avail Portal	20
Multiparty Events   Setup	
Multiparty Events   Joining	



# **Remote User | Initial Setup**

A Remote User has the ability to:

- ✓ Schedule availability
- ✓ Accept incoming calls
- Update user profile
- ✓ Request facility access

Remote User Admins have the added ability to:

- ✓ Invite new users
- ✓ View/edit team's availability

### **Initial Remote App Setup**

#### **Download Remote App**

- Search "Avail Medsystems" in App Store on your iPad
- Select "Avail Remote App"; tap "Get"
- Launch App
  - Read & accept EULA
- Login using the same credentials as the Avail Portal
  - Select notification preferences when prompted
  - At 1st login, grant app camera & microphone access





# **Remote User | Avail Portal**

### Scheduler



### **Facility Request Status**

User profile	General information		Request a to a facilit	access y •·····	Add new Facility	
Nadia Hakeem	Therapies & procedures	Name	State	Request status	Actions	
Manager	Preferences	WellState Hospital	Pennsylvania	Approved	:	
Select Manager 💌		Sea Cove Hospital	New York	Approved	:	
		River Oaks Hospital	Arkansas	Approved	:	
				$\bigcirc$		
				Ť		
				<b>i</b>		
			View re	equest		
	Nadia Hakeem Manager Select Mänager	Nadia Hakeem   Manager   Select Manager   Preferences	Nadia Hakeem     Therapies & procedures       Manager     Facilities     Name       Select Manager     Preferences     WellState Hospital       Sea Cove Hospital     River Oaks Hospital	Nadia Hakeem       Therapies & procedures         Manager       Facilities       Name       State         Select Manager       Preferences       WellState Hospital       Pennsylvania         Sea Cove Hospital       New York       River Oaks Hospital       Arkansas         View red       View red       View red	Nadia Hakeem       Therapies & procedures         Manager       Facilities       Name       State       Request status         Select Manager       Preferences       WellState Hospital       Pennsylvania       Approved         Select Manager       River Oaks Hospital       New York       Approved         Well State Hospital       New York       Approved         Kiver Oaks Hospital       Arkansas       Approved         View request       View request	



# **Remote User | Avail Portal**

# **Incoming Facility Call**

🖌 Avail			Nadia Hakeem [	Vendor Field Representa	tive, Peripheral Vascular)
SCHEDULER	NH Nadia Hakeem			Timezone	America/Los_Angeles
	Today < > July 12 - 18			n	nonth week day
	12 Answer incoming			FRI 17	sat 18
	call from Facility 7:30am - 6:00p	Sea Cove Hospital is calling you.	m	7:30am - 6:00pm	
	9am	ANSWER DECLINE			
	10am				
	llam				
	12pm				

# **Missed Facility Call**

10:55 AM	Spine, Cervical Fusions   Northstar Medical Center, Quincy MA	CALL	васк	🖊 Av	vail	Symmetry Surgical Devices   Marcus Thompson
🖌 Avail		Symmetry Surgical Devices   Marcus Thor	rpson M	•	Back > Missed Calls	
•	Marcus Thompson	Timezone: Los Angele Part Misse	I Calls		MT Marcus Thompson	
	DDAY C > January 26 - February 01	MONTH OTD My O	panization		TODAY	Temezone: Los Angeles - Pacific 🗸
н	UR SUN 26 MON 27 TUE 28 WED	29 THU 30 FRI 31 SAT	01		WED 29	
1						
2	m				10:55 AM	CALL BACK
3					Hybrid 1	
-					Quincy, MA	
6		View missed				
7	m				9:27 AM EP. Ablation	
8	2:00 am - 11:30 am	Calls			Cath Lab 7 Balla Virta Surgical Contor	A blue "CALL BACK" button
9	8.30 am - 100 pm				Mount Hamilton, CT	indicates a call back to the
10	930 am - 200 pm	9:30 am- 2.00 pm				
	m				9:05 AM	Console may be possible
12					Cath Lab 7 Rolls Vista Surgical Conter	
					Mount Hamilton, CT	
3					TUF 28	
4	100 pm - 100 pm					
5	an a				5:40 PM	
6	pas				Spine, Lateral Fusions Lab 2	
7	m				Sanders Medical Bangor, ME	
	ha la					
9					2:32 PM	CALL BACK
10					Coronary, Heart Catheterization Cathether Lab	
					Chamberlain Health South Portland, ME	



# **Multiparty Events | Setup**

### **Creating an Event in the Avail Portal**

Avail						Jack	Black [HCP, Demo-hcp]
SCHEDULER	JB Jack Black Today <	Click "Events" to create or see a list of events	тие 20	wed 21	тни 22	Timezone m	America/Los_Angeles * onth week day SAT 24
	7am						
	Bam						

clive Events					
oril 23		No Events to	join	Click "Create Event"; ensure you are on the today's date	CREATE EVEN
Event Title		1			
TAVI in Low Risk Patient with Aortic Stenosis			C	omplete	
Facility			e	vent	
WellSpan Surgeon/Specialist	-		in	formation	
Pedro Blanca, MD			a	nd click	
Event Host sruthi hcp (sruthi-hcp@qa.com)	÷		" (	Submit"	
Attendee 1	_				
petesampas@kingcardiovascularllc.com					
Attendee 2					
alisejohnson@kingcardiovascularllc.com					
Attendee 3					
jane.dempsey@vlchealthsystem.org					



# **Multiparty Events | Joining**

### **Joining via Avail Portal**

🖊 Avail					Jack 8	Black [HCP, Demo-hcp]
SCHEDULER					Timezone	America/Los_Angeles *
Today <	Click "Events" to see a list of events	тие 20	WED 21	тни 22	FRI 23	sat 24
< Return to Portal Active Events April 23 Pedro Blanca, MD TAVI in Low Risk Patient v WellSpan Jack Black	Locate event			Click "Join	n" •	

# Joining via Email Notification

### Starting/Joining via Avail Console

	Avail 1 0
	CATH LAB 2 EVENTS JANUARY 18, 2021
Avail	Events  C. Boh Hame EndoVive Enteral Access Initial Placement: Standard Perculaneous Endoscepic Gastrostomy (PEG) Kit Containeous
The future of collaboration is here.	EP Coronary Stoch Jameson TAVR Minimally Invasive Alternative for Storage tic Stoch Jameson TAVR Minimally Invasive Alternative for Storage tic
[Tom Phan] is inviting you to join an Avail Event Case: DODOX(] When: DODOX]	From the
Broadcasting from: [WellSpan - Adams] Surgeon: [Dr. X] Hosted by: [John Dow]	A page, select
Click "Join	Events Event to Start
(http://event)	Surgeon Operator Event Title and Status Host Name